



QUICK START GUIDE

Learn how to use and store
once-weekly Ozempic®

This material is a **DRAFT** for preparatory **internal use in NN only**. Affiliates are responsible for reviewing the promotional material against local label, more stringent relevant local legislation, and if relevant local code of conduct (cf. S.O.P 100965 Approval of Promotional Material in Novo Nordisk) before distribution.

Please see the package leaflet for more information or contact your doctor or nurse with any questions about how to manage your type 2 diabetes or use Ozempic®.



ONCE-WEEKLY
OZEMPIC®
semaglutide injection

How do I use the once-weekly Ozempic® pen?

A Pen cap

Always replace it to protect Ozempic® from light

B Pen window

Look through it to check that Ozempic® is clear and colourless

C Dose counter

Use it to see the dose you select

D Dose selector

Turn it to select your dose

E Dose button

Press and hold it to inject your dose



Where to inject

The best places on your body to inject Ozempic® are the front of your thighs, the front of your waist (abdomen) or your upper arms.


How much to inject

The starting dose of Ozempic® is 0.25 mg for the first 4 weeks of treatment.

How to store

After opening, Ozempic® may be stored for 6 weeks at a temperature below 30°C or in a refrigerator at 2°C to 8°C.

Please see the package leaflet for more information or contact your doctor or nurse with any questions about how to manage your type 2 diabetes or use Ozempic®.



Injection steps

Please see the full instructions for use that came with your Ozempic® pen.



1 Check your pen and attach a new needle

Make sure your pen contains **clear and colourless** Ozempic®. Screw on a new needle and pull off **both** needle caps.



2 Check the flow

With each **new** pen, turn the dose selector to select the **flow check symbol** (▬ ▬ ▬). Press and hold in the dose button. Repeat until a drop appears.




3 Select your dose

Turn the dose selector until the dose counter shows your dose (0.25 mg, 0.5 mg or 1 mg).



4 Inject your dose

Insert the needle. Press and hold down the dose button. After the dose counter reaches 0, **slowly count to 6**.



5 Remove the needle

Carefully remove and dispose of the needle. Put the pen cap back on.

This material is a **DRAFT** for preparatory **internal use in NN only**. Affiliates are responsible for reviewing the promotional material against local label, more stringent relevant local legislation, and if relevant local code of conduct (cf. S.O.P 100965 Approval of Promotional Material in Novo Nordisk) before distribution.

This material is a **DRAFT** for preparatory **internal use in NN only**. Affiliates are responsible for reviewing the promotional material against local label, more stringent relevant local legislation, and if relevant local code of conduct (cf. S.O.P 100965 Approval of Promotional Material in Novo Nordisk) before distribution.



Please see the package leaflet for more information or contact your doctor or nurse with any questions about how to manage your type 2 diabetes or use Ozempic®.

Ozempic® is a registered trademark of Novo Nordisk A/S.
2018 © Novo Nordisk A/S, Novo Allé, DK-2880, Bagsværd, Denmark
HQMM/SM/1117/0132(1), Approval date: October 2018

ONCE-WEEKLY
OZEMPIC®
semaglutide injection